



Making a **World of Difference,**
to a **World of People.**



Product Catalog

Unistik® 3


Side-Activated Safety Lancets

- **Greater Comfort**
Features Comfort Zone Technology and side-activation, designed to reduce pain.¹
- **Safely Protects**
Permanently retracts to help protect against accidental needlestick and cross contamination.
- **Exceptional Quality**
Made in England to exacting quality standards.






Clinical Evidence

Lancet trials showed that **Comfort Zone Technology® (CZT®)** reduced the overall pain experienced during the blood sampling procedure when compared to devices without **Comfort Zone Technology®**.¹

Maximum Pain Score comparison with and without CZT®



Clinical Evaluation results, March 2007

Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Unistik® 3 Gentle 30G, 1.5mm For frequent testing 	AT 1028	08470-1028-01	3-8470102801-4	25	100	3.5" x 3.25" x 2.375"	0.2 lbs	17" x 12.5" x 14"	21.1 lbs
	AT 1027	08470-1027-01	3-8470102701-7	50	100	4" x 2.75" x 3.5"	0.4 lbs	17.25" x 14.25" x 19.25"	40 lbs
	AT 1022	08470-1022-01	3-8470102201-2	100	40	4" x 3.625" x 5.25"	0.7 lbs	20.5" x 11.25" x 15.25"	31 lbs
	AT 1024	08470-1024-01	3-8470102401-6	200	20	5.5" x 4.5" x 6"	1.4 lbs	22.75" x 12.25" x 11.5"	29 lbs
Unistik® 3 Comfort 28G, 1.8mm Low Blood Flow 	AT 1048	08470-1048-01	3-8470104801-2	25	100	3.5" x 3.25" x 2.375"	0.2 lbs	17" x 12.5" x 14"	21.1 lbs
	AT 1047	08470-1047-01	3-8470104701-5	50	100	4" x 2.75" x 3.5"	0.4 lbs	17.25" x 14.25" x 19.25"	40 lbs
	AT 1042	08470-1042-01	3-8470104201-0	100	40	4" x 3.625" x 5.25"	0.7 lbs	20.5" x 11.25" x 15.25"	31 lbs
	AT 1044	08470-1044-01	3-8470104401-4	200	20	5.5" x 4.5" x 6"	1.4 lbs	22.75" x 12.25" x 11.5"	29 lbs
Unistik® 3 Normal 23G, 1.8mm Medium Blood Flow 	AT 1008	08470-1008-01	3-8470100801-6	25	100	3.5" x 3.25" x 2.375"	0.2 lbs	17" x 12.5" x 14"	21.1 lbs
	AT 1007	08470-1007-01	3-8470100701-9	50	100	4" x 2.75" x 3.5"	0.4 lbs	17.25" x 14.25" x 19.25"	40 lbs
	AT 1002	08470-1002-01	3-8470100201-4	100	40	4" x 3.625" x 5.25"	0.7 lbs	20.5" x 11.25" x 15.25"	31 lbs
	AT 1004	08470-1004-01	3-8470100401-8	200	20	5.5" x 4.5" x 6"	1.4 lbs	22.75" x 12.25" x 11.5"	29 lbs
Unistik® 3 Extra 21G, 2.0mm High Blood Flow 	AT 1018	08470-1018-01	3-8470101801-5	25	100	3.5" x 3.25" x 2.375"	0.2 lbs	17" x 12.5" x 14"	21.1 lbs
	AT 1017	08470-1017-01	3-8470101701-8	50	100	4" x 2.75" x 3.5"	0.4 lbs	17.25" x 14.25" x 19.25"	40 lbs
	AT 1012	08470-1012-01	3-8470101201-3	100	40	4" x 3.625" x 5.25"	0.7 lbs	20.5" x 11.25" x 15.25"	31 lbs
	AT 1014	08470-1014-01	3-8470101401-7	200	20	5.5" x 4.5" x 6"	1.4 lbs	22.75" x 12.25" x 11.5"	29 lbs
Unistik® 3 Dual 18G, 1.8mm Maximum flow fingersticks and low flow heelsticks 	AT 1062	08470-1062-01	3-8470106201-8	100	40	4" x 3.625" x 5.25"	0.7 lbs	20.5" x 11.25" x 15.25"	31 lbs
	AT 1064	08470-1064-01	3-8470106401-2	200	20	5.5" x 4.5" x 6"	1.4 lbs	22.75" x 12.25" x 11.5"	29 lbs

¹Data on file.

Unistik® Touch™


Contact-Activated Safety Lancets

- Enhanced Comfort**
Features Comfort Zone Technology, designed to reduce pain.¹
- Safely Protects**
Permanently retracts to help protect against accidental needlestick and cross contamination.
- Easy-to-Use**
Intuitive design requires minimal training and can be used in three easy steps.

Clinical Evidence

Lancet trials showed that **Comfort Zone Technology® (CZT®)** reduced the overall pain experienced during the blood sampling procedure when compared to devices without **Comfort Zone Technology®**.¹

Maximum Pain Score comparison with and without CZT®



Clinical Evaluation results, March 2007

Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Unistik® Touch™ 30G, 1.5mm Lowest Blood Flow	AT 1422	08470-1422-01	3-8470142201-0	100	40	4.875" x 2.625" x 3.25"	0.45 lbs	17.25" x 9.75" x 11.5"	19.65 lbs
	AT 1424	08470-1424-01	3-8470142401-4	200	20	4" x 3.25" x 5.25"	0.9 lbs	17.25" x 9.75" x 11.5"	19.65 lbs
Unistik® Touch™ 28G, 1.8mm Low Blood Flow	AT 1442	08470-1442-01	3-8470144201-8	100	40	4.875" x 2.625" x 3.25"	0.45 lbs	17.25" x 9.75" x 11.5"	19.65 lbs
	AT 1444	08470-1444-01	3-8470144401-2	200	20	4" x 3.25" x 5.25"	0.9 lbs	17.25" x 9.75" x 11.5"	19.65 lbs
Unistik® Touch™ 23G, 2.0mm Medium Blood Flow	AT 1402	08470-1402-01	3-8470140201-2	100	40	4.875" x 2.625" x 3.25"	0.45 lbs	17.25" x 9.75" x 11.5"	19.65 lbs
	AT 1404	08470-1404-01	3-8470140401-6	200	20	4" x 3.25" x 5.25"	0.9 lbs	17.25" x 9.75" x 11.5"	19.65 lbs
Unistik® Touch™ 21G, 2.0mm High Blood Flow	AT 1412	08470-1412-01	3-8470141201-1	100	40	4.875" x 2.625" x 3.25"	0.45 lbs	17.25" x 9.75" x 11.5"	19.65 lbs
	AT 1414	08470-1414-01	3-8470141401-5	200	20	4" x 3.25" x 5.25"	0.9 lbs	17.25" x 9.75" x 11.5"	19.65 lbs



¹Data on file.

Unistik® Pro

Top-Activated Safety Lancets

- **Enhanced Comfort**

Features Comfort Zone Technology, designed to reduce pain.¹




- **Safely Protects**

Permanently retracts to help protect against accidental needlestick and cross contamination.

- **Easy-to-Use**

Featuring a syringe-like grip and a simple push-button design, Unistik Pro makes blood sampling easy and intuitive for healthcare professionals.



Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Unistik® Pro 28G, 1.2mm Low Blood Flow 	AT 1644	08470-1644-01	3-8470164401-6	200	20	5.5" x 3.5" x 4.75"	0.75 lbs	11.42" x 24" x 7.7"	15 lbs
Unistik® Pro 25G, 1.6mm Medium Blood Flow 	AT 1634	08470-1634-01	3-8470163401-7	200	20	5.5" x 3.5" x 4.75"	0.75 lbs	11.42" x 24" x 7.7"	15 lbs
Unistik® Pro 21G, 2.0mm High Blood Flow 	AT 1614	08470-1614-01	3-8470161401-9	200	20	5.5" x 3.5" x 4.75"	0.75 lbs	11.42" x 24" x 7.7"	15 lbs






Unilet® ComforTouch®

Home Use Recappable Lancets

- Recappable design for safer disposal.
- Precision-sharpened for enhanced comfort.
- Universal fit with most lancing devices.¹



Unilet® Lancets

Home Use Lancets

- Cost effective solution for a low cash price.
- Precision-sharpened for enhanced comfort.
- Universal fit with most lancing devices.¹

ComforTouch 30G

ComforTouch 28G

ComforTouch 26G

ComforTouch 23G



Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Unilet® ComforTouch® 30G Super Thin Lancets	AT 0465	08470-0465-01	3-8470046501-8	100	200	2.625" x 3.25" x 1.625"	0.1 lbs	17.25" x 12.75" x 14"	24.3 lbs
	AT 0460	08470-0460-01	3-8470046001-3	200	100	3.375" x 3.25" x 2.375"	0.25 lbs	17.25" x 12.75" x 14"	23.4 lbs
Unilet® ComforTouch® 28G Ultra Thin Lancets	AT 0925	08470-0925-01	3-8470092501-7	100	200	2.625" x 3.25" x 1.625"	0.1 lbs	17.25" x 12.75" x 14"	24.3 lbs
	AT 0920	08470-0920-01	3-8470092001-2	200	100	3.375" x 3.25" x 2.375"	0.25 lbs	17.25" x 12.75" x 14"	23.4 lbs
Unilet® ComforTouch® 26G Alternate Site Lancets	AT 0435	08470-0435-01	3-8470043501-1	100	200	2.625" x 3.25" x 1.625"	0.1 lbs	17.25" x 12.75" x 14"	24.3 lbs
	AT 0430	08470-0430-01	3-8470043001-6	200	100	3.375" x 3.25" x 2.375"	0.25 lbs	17.25" x 12.75" x 14"	23.4 lbs
Unilet® ComforTouch® 23G Thin Lancets	AT 0455	08470-0455-01	3-8470045501-9	100	200	2.625" x 3.25" x 1.625"	0.1 lbs	17.25" x 12.75" x 14"	24.3 lbs
	AT 0450	08470-0450-01	3-8470045001-4	200	100	3.375" x 3.25" x 2.375"	0.25 lbs	17.25" x 12.75" x 14"	23.4 lbs

Unilet 33G

Unilet 30G

Unilet 28G



Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Unilet® 33G Micro Thin Lancets	AT 0585	08470-0585-01	3-8470058501-3	100	200	2.625" x 3.25" x 1.625"	0.1 lbs	17.25" x 12.75" x 14"	24.3 lbs
Unilet® 30G Super Thin Lancets	AT 0575	08470-0575-01	3-8470057501-4	100	200	2.625" x 3.25" x 1.625"	0.1 lbs	17.25" x 12.75" x 14"	24.3 lbs
Unilet® 28G Ultra Thin Lancets	AT 0565	08470-0565-01	3-8470056501-5	100	200	2.625" x 3.25" x 1.625"	0.1 lbs	17.25" x 12.75" x 14"	24.3 lbs


¹Not compatible with ACCU-CHEK® Softclix or OneTouch® Delica™



Autolet®

Advanced Lancing Device


- Automatically primes the lancet upon insertion.
- Easy cap removal and lancet release.
- Features Comfort Zone Technology designed for maximum patient comfort¹.
- 9 adjustable settings.



Clinical Evidence

Lancet trials showed that **Comfort Zone Technology® (CZT®)** reduced the overall pain experienced during the blood sampling procedure when compared to devices without **Comfort Zone Technology®**.¹

Maximum Pain Score comparison with and without CZT®



Clinical Evaluation results, March 2007

Autolet®

Lancing Device




Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Autolet® Lancing Device (Box)	AT 0270	08470-0270-01	3-8470027001-8	1	300	2.125" x 1.375" x 3.625"	0.05 lbs	11.75" x 15.5" x 21"	18.8 lbs
Autolet® Lancing Device (Bag)	AT 0271	08470-0271-01	3-8470027101-5	1	100	Bag	0.05 lbs	12.5" x 7.5" x 17.25"	6.25 lbs

Unistik® TinyTouch™

Heel Incision Device

- Designed for Comfort**
Sweeping motion blade designed to help reduce discomfort during sampling.
- Easy to Use**
Contoured, ergonomic design provides versatile grips with alignment guides for precise placement.
- Safely Protects**
Permanently retracts to help protect against accidental incisions and cross contamination.

Incision comparison in silicone¹



Unistik® TinyTouch™ vs Leading Competitor



Premie Heel Incision Device

0.85mm (D) x 1.75mm (L)



Full Term Heel Incision Device

1.00mm (D) x 2.50mm (L)



18G Single Use Safety Lancet

1.80mm (D) x 18G

Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Unistik® TinyTouch® Full Term Depth 1.0mm x Length 2.5mm	AT 1250	08470-1250-01	3-8470125001-9	50	30	4" x 3.56" x 4.38"	0.55 lbs	22.25" x 8.875" x 13.875"	19.45 lbs
Unistik® TinyTouch® Premie Depth .85mm x Length 1.75mm	AT 1200	08470-1200-01	3-8470120001-4	50	30	4" x 3.56" x 4.38"	0.55 lbs	22.25" x 8.875" x 13.875"	19.45 lbs
Unistik® TinyTouch® Single Use Safety Lancet Depth 1.8mm x 18G	AT 1267	08470-1267-01	3-8470126701-7	50	100	4" x 2.75" x 3.5"	0.4 lbs	17.25" x 14.25" x 19.25"	40 lbs
Unistik® TinyTouch® Full Term: Individually Bagged Depth 1.0mm x Length 2.5mm	AT 1260	08470-1260-01	3-8470126001-8	25	24	4" x 4.375" x 3.625"	0.4 lbs	22.25" x 8.875" x 13.875"	11.85 lbs
Unistik® TinyTouch® Premie: Individually Bagged Depth .85mm x Length 1.75mm	AT 1210	08470-1210-01	3-8470121001-3	25	24	4" x 4.375" x 3.625"	0.4 lbs	22.25" x 8.875" x 13.875"	11.85 lbs

¹Image depicted shows magnification of incisions made in silicone comparing TinyTouch versus a leading national brand. (2013).



Autoject® 2

Automatic Injection Aid

- **Hides needle to reduce anxiety**
Designed to conceal the needle from the user's view prior to injecting, reducing anxiety for many patients who are frightened by the sight of needles.
- **Easy to use**
Ergonomic design allows patients to self-inject, reducing their reliance on others. Especially beneficial for children and seniors.
- **Safety interlock mechanism**
Designed to help prevent accidental firing until the device is pressed against the skin.

Autoject® 2 is not for use with glass syringes.

Fixed needle



Removable needle



EI Fixed Needle



Product Description	Item Code	NDC	UPC	Box Qty	Boxes/ Case	Item Size	Item Weight	Case Size	Case Weight
Autoject® 2 Fixed needle	AJ 1300	08470-1300-01	3-8470130001-1	1	50	4.875" x 6.75" x 1.5"	0.35 lbs	20.75" x 17" x 10.5"	17.75 lbs
Autoject® 2 EI Fixed Needle	AJ 1310	08470-1310-01	3-8470131001-0	1	50	4.875" x 6.75" x 1.5"	0.35 lbs	20.75" x 17" x 10.5"	17.75 lbs
Autoject® 2 Removable needle	AJ 1311	08470-1311-01	3-8470131101-7	1	50	4.875" x 6.75" x 1.5"	0.35 lbs	20.75" x 17" x 10.5"	17.75 lbs



Autoject® EI has a built-in indicator to help individuals with visual impairments know when their dose is complete.



Compatible Syringes

The Autoject 2 may be used for subcutaneous injection with the following syringes:

Syringe Brand	Size	Needle Length	Depth Adjuster	Depth Adjuster Setting (6mm/ 1/4")
B. Braun	1.0ml (1cc)	8mm (5/16")	Long (blue scale)	c
B. Braun	0.5ml (1/2cc)	8mm (5/16")	Long (blue scale)	n
B. Braun	1.0ml (1cc)	12mm (1/2")	Long (blue scale)	e
B. Braun	0.5ml (1/2cc)	12mm (1/2")	Long (blue scale)	q
BD Micro-Fine®/Ultra-Fine®	1.0ml (1cc)	12.7mm (1/2")	Long (blue scale)	c
BD Micro-Fine®/Ultra-Fine®	0.5ml (1/2cc)	8mm (5/16")	Short (black scale)	m
BD Micro-Fine®/Ultra-Fine®	0.5ml (1/2cc)	12.7mm (1/2")	Short (black scale)	o
BD Micro-Fine®/Ultra-Fine®	0.3ml (3/10cc)	12.7mm (1/2")	Short (black scale)	j
Precision Sure-Dose™/Terumo®	1.0ml (1cc)	12mm (1/2")	Long (blue scale)	b
Precision Sure-Dose™/Terumo®	0.5ml (1/2cc)	12mm (1/2")	Short (black scale)	l
Precision Sure-Dose™/Terumo®	0.3ml (3/10cc)	12mm (1/2")	Short (black scale)	l

Brands are trademarks of their respective owners.

Autopen®

Insulin Drug Delivery Pens

- **Simple, automatic delivery**
Using a unique side-mounted release button.
- **Easy to use**
Dose selector with audible clicks.
- **Compatible with most pen needles¹**
Including Unifine® Pentips® and Unifine® Pentips® Plus

Accessories included with Autopen®



Release Button Extension
Makes it easier to push the release button.



Dose Selector Adaptor
Makes dialing a dose easier.

Autopen® Classic




1-21 units



2-42 units

Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Autopen® Classic 1-21 units in increments of 1 unit Compatible with Lilly insulin	AN 3810	08470-3810-01	3-8470381001-3	1	100	3.25" x 7.75" x 1.25"	0.25 lbs	22.5" x 22.75" x 10"	27.45 lbs
Autopen® Classic 2-42 units in increments of 2 unit Compatible with Lilly insulin	AN 3800	08470-3800-01	3-8470380001-4	1	100	3.25" x 7.75" x 1.25"	0.25 lbs	22.5" x 22.75" x 10"	27.45 lbs

¹Data on file.



Unifine® Pentips® Plus

Pen Needles with a Built-In Remover

- **The first and only pen needle with a built-in remover**
Designed to be easier, safer and more convenient¹.
- **Ergonomic design**
Designed to make pen needle attachment and removal easy, with a larger container for improved grip and control.
- **Comprehensive coverage**
Covered under most insurance plans.²

Drug Delivery



Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Unifine® Pentips® Plus 0.23mm (33G) x 4mm	AN 3860	08470-3860-01	3-8470386001-8	100	24	4" x 5.25" x 3.625"	0.45 lbs	17.5" x 12" x 12.75"	12.85 lbs
	AN 1860	08470-1860-01	3-8470186001-0	30	100	3.375" x 3.25" x 2.375"	0.1 lbs	17.25" x 12.75" x 14"	20.5 lbs
Unifine® Pentips® Plus 0.23mm (32G) x 4mm	AN 3840	08470-3840-01	3-8470384001-0	100	24	4" x 5.25" x 3.625"	0.45 lbs	17.5" x 12" x 12.75"	12.85 lbs
	AN 1840	08470-1840-01	3-8470184001-2	30	100	3.375" x 3.25" x 2.375"	0.1 lbs	17.25" x 12.75" x 14"	20.5 lbs
Unifine® Pentips® Plus 0.25mm (31G) x 5mm	AN 3850	08470-3850-01	3-8470385001-9	100	24	4" x 5.25" x 3.625"	0.45 lbs	17.5" x 12" x 12.75"	12.85 lbs
	AN 1850	08470-1850-01	3-8470185001-1	30	100	3.375" x 3.25" x 2.375"	0.1 lbs	17.25" x 12.75" x 14"	20.5 lbs
Unifine® Pentips® Plus 0.25mm (31G) x 6mm	AN 3890	08470-3890-01	3-8470389001-5	100	24	4" x 5.25" x 3.625"	0.45 lbs	17.5" x 12" x 12.75"	12.85 lbs
	AN 1890	08470-1890-01	3-8470189001-7	30	100	3.375" x 3.25" x 2.375"	0.1 lbs	17.25" x 12.75" x 14"	20.5 lbs
Unifine® Pentips® Plus 0.25mm (31G) x 8mm	AN 3830	08470-3830-01	3-8470383001-1	100	24	4" x 5.25" x 3.625"	0.45 lbs	17.5" x 12" x 12.75"	12.85 lbs
	AN 1830	08470-1830-01	3-8470183001-3	30	100	3.375" x 3.25" x 2.375"	0.1 lbs	17.25" x 12.75" x 14"	20.5 lbs
Unifine® Pentips® Plus 0.33mm (29G) x 12mm	AN 3829	08470-3829-01	3-8470382901-5	100	24	4" x 5.25" x 3.625"	0.45 lbs	17.5" x 12" x 12.75"	12.85 lbs
	AN 1829	08470-1829-01	3-8470182901-7	30	100	3.375" x 3.25" x 2.375"	0.1 lbs	17.25" x 12.75" x 14"	20.5 lbs

Also available in 90ct.



Also Available



¹HRW (2014) Impact of Unifine Pentips Plus on pen needle changing behaviour amongst people with diabetes medicating with injectable formats.

²Independent and chain pharmacy co-pay adjudication (April 2016). Actual coverage and co-pay may vary from setting to setting, and insurer to insurer. Data on file.

Unifine® Pentips®

Pen Needles

- **Designed with patient comfort in mind**
Featuring OptiFlow and DiamondPoint Technology for low glide and penetration force¹.
- **Universal fit**
Compatible with most injection pens².
- **Comprehensive coverage**
Covered under most insurance plans.³



Also Available



**Unifine® Pentips®
Remover**

Available for individual sale.

Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Unifine® Pentips® 0.23mm (33G) x 4mm	AN 3560	08470-3560-01	3-8470356001-7	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs
	AN 1160	08470-1160-01	3-8470116001-1	30	50	2.625" x 3.25" x 1.625"	0.1 lbs	12" x 9" x 7"	5.1 lbs
Unifine® Pentips® 0.23mm (32G) x 4mm	AN 3540	08470-3540-01	3-8470354001-9	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs
	AN 1140	08470-1140-01	3-8470114001-3	30	50	2.625" x 3.25" x 1.625"	0.1 lbs	12" x 9" x 7"	5.1 lbs
Unifine® Pentips® 0.25mm (31G) x 5mm	AN 3550	08470-3550-01	3-8470355001-8	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs
	AN 1150	08470-1150-01	3-8470115001-2	30	50	2.625" x 3.25" x 1.625"	0.1 lbs	12" x 9" x 7"	5.1 lbs
Unifine® Pentips® 0.25mm (31G) x 6mm	AN 3590	08470-3590-01	3-8470359001-4	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs
	AN 1190	08470-1190-01	3-8470119001-8	30	50	2.625" x 3.25" x 1.625"	0.1 lbs	12" x 9" x 7"	5.1 lbs
Unifine® Pentips® 0.25mm (31G) x 8mm	AN 3530	08470-3530-01	3-8470353001-0	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs
	AN 1130	08470-1130-01	3-8470113001-4	30	50	2.625" x 3.25" x 1.625"	0.1 lbs	12" x 9" x 7"	5.1 lbs
Unifine® Pentips® 0.33mm (29G) x 12mm	AN 3529	08470-3529-01	3-8470352901-4	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs

Also available in 90ct.





PENTIPS®

Pen Needles

- Quality pen needles at an affordable price.
- Universal fit with most injection pens.¹
- Covered under most insurance plans.²

Product Description	Item Code	NDC	UPC	Box Qty	Boxes/ Case	Item Size	Item Weight	Case Size	Case Weight
Pentips® 0.23mm (32G) x 4mm	AN 3440	08470-3440-01	3-8470344001-2	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs
Pentips® 0.25mm (31G) x 5mm	AN 3450	08470-3450-01	3-8470345001-1	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs
Pentips® 0.25mm (31G) x 6mm	AN 3490	08470-3490-01	3-8470349001-7	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs
Pentips® 0.25mm (31G) x 8mm	AN 3430	08470-3430-01	3-8470343001-3	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs
Pentips® 0.33mm (29G) x 12mm	AN 3429	08470-3429-01	3-8470342901-7	100	50	4" x 2.75" x 3.5"	0.3 lbs	21.75" x 14.25" x 9.375"	17.6 lbs

¹Data on file. ²Actual coverage and co-pay may vary from setting to setting, and insurer to insurer. Data on file.



AutoDrop®

Eye Drop Guide

- Perfect solution for conserving expensive eye drops.
- Fits most eye drop bottles on the market.
- Holds eye open and directs medication to help reduce waste.

AutoSqueeze™

Bottle Aid

- Makes bottle squeezing easier.
- Clips on eye drop bottle to provide extra leverage.
- Excellent solution for people with arthritis or dexterity issues.

AutoDrop®



AutoSqueeze®



Product Description	Item Code	NDC	UPC	Box Qty	Boxes/ Case	Item Size	Item Weight	Case Size	Case Weight
AutoDrop® Eye Drop Guide	OP 6000	08470-6000-01	3-8470600001-5	1	200	2.25" x 1.375" x 4.25"	0.05 lbs	19.5" x 12.5" x 15.75"	12.5 lbs
AutoSqueeze® Bottle Aid	OP 6100	08470-6100-01	3-8470610001-2	1	200	2.25" x 1.375" x 4.25"	0.05 lbs	19.5" x 12.5" x 15.75"	12.5 lbs



Neuropen[®]

Neuropathy Screening Tools

- **Effective Screening Aid**

Combines two calibrated tests, enabling clinicians to identify patients at risk of foot ulceration.

1. **10g Pressure Test**

The Neuropen Monofilaments provide a measured 10g force designed to help assess reduced pressure perception.

2. **40g Sharpness Test**

The Neurotip provides a measured 40g force designed to help assess reduced sensation.

- **Replaceable**


The Neuropen Monofilaments and Neurotip can be replaced as required.

Neuropen[®]

Monofilaments

Monofilaments test for assessing protective touch/pressure sensation in large nerve fibers.

Suggested Neuropathy Testing Sites



Left foot

Right foot

Neuropen[®]



Neurotips[™]

Neurotips test for assessing reduced sensation to sharpness/pain in small nerve fibers. The calibrated spring mechanism exerts a force of 40g to help assess patients.



Product Description	Item Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Neuropen[®] Neuropathy Screening Device	NT 0100	08470-0100-01	3-8470010001-8	1	100	7.25" x 2.75" x 1"	0.05 lbs	12.5" x 7.5" x 17.25"	6.65 lbs
Neuropen[®] Monofilaments Replacement Monofilaments	NT 0104	08470-0104-01	3-8470010401-6	5	100	3.75" x 3.5" x 1"	0.1 lbs	16.875" x 13.25" x 11.25"	14.15 lbs
Neurotips[™] Neurological Examination Pins	NT 5405	08470-5405-01	3-8470540501-9	100	100	3" x 1.875" x 2.625"	0.15 lbs	16.875" x 13.25" x 11.25"	20 lbs

Rapport™

Erectile Dysfunction Management

- Vacuum therapy system helps to quickly produce an erection by drawing blood into the penis.
- Designed to be a safer, reliable alternative to prescription ED medication.¹
- Long-term use of this device can equate to significant cost savings versus expensive prescription drug options.
- Easy to use with up to seven ring options.
- Compact kit for storage or to take on-the-go.



The Rapport™ Classic system includes:

- 1 clear plastic vacuum tube
- 1 hand pump with a pressure release valve
- 5 sizes of penile rings (sizes 4-8)
- 2 sizes of transfer sleeves
- 1 loading cone
- 1 tube of water-soluble lubricating jelly
- 1 demonstration DVD
- Patient instructions for use
- All the above items are supplied in a discreet carrying case.

Product Description	Code	NDC	UPC	HCPCS Code	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Rapport™ Classic	SM 2000	08470-2000-01	3-8470200001-9	L7902	1	5	10.625" x 5" x 6.875"	1.85 lbs	17.25" x 12.75" x 14"	10 lbs
Rapport™ Rings (Size 3-9)	Size 3 SM 2203	08470-2203-01	3-8470220301-4	L7902	1	–	3" x 5" x 0.25"	0.05 lbs	–	–
	Size 4 SM 2204	08470-2204-01	3-8470220401-1							
	Size 5 SM 2205	08470-2205-01	3-8470220501-8							
	Size 6 SM 2206	08470-2206-01	3-8470220601-5							
	Size 7 SM 2207	08470-2207-01	3-8470220701-2							
	Size 8 SM 2208	08470-2208-01	3-8470220801-9							
	Size 9 SM 2209	08470-2209-01	3-8470220901-6							
Rapport™ Pump Handle	SM 2009	08470-2009-01	3-8470200901-2	L7902	1	10	6" x 8" x 2.5"	0.20 lbs	11" x 8" x 6"	2.35 lbs

¹International Journal of Impotence Research (2003) 15, 63–71.

Sexual Wellness Solutions for Women



Size 0 - 7cm x 1.5cm (2.75" x 0.6")



Size 1 - 9cm x 2cm (3.54" x 0.8")




Size 2 - 11cm x 2.5cm (4.33" x 1")



Size 3 - 14cm x 3cm (5.51" x 1.2")



Size 4 - 16cm x 3.5cm (6.30" x 1.4")

Twist-and-lock
universal handle

Amielle® Restore

Vaginal Exercise System

Product Description		Item Code	NDC	UPC	HCPCS Code	Items Qty	Items/ Case	Item Size	Item Weight	Case Size	Case Weight
Amielle® Comfort Rx-Only Vaginal Dilator System Gentle self treatment for vaginismus and dyspareunia		SM 2100	08470-2100-01	3-8470210001-6	E1399	1	25	7.25" x 4.25" x 3"	0.6 lbs	19.5" x 12.5" x 15.75"	15.9 lbs
Amielle® Restore OTC Vaginal Exercise System Designed to address pain and discomfort during penetration		SM 3100	08470-3100-01	3-8470310001-5	E1399	1	25	7.25" x 4.25" x 3"	0.6 lbs	19.5" x 12.5" x 15.75"	15.9 lbs
Amielle® Dilators Rx Only (Size 0-4)	Size 0	SM 2101	08470-2101-01	3-8470210101-3	E1399	1	50	6" x 8" x 1"	0.1 lbs	21.75" x 14.25" x 9.375"	5.5 lbs
	Size 1	SM 2102	08470-2102-01	3-8470210201-0							6.2 lbs
	Size 2	SM 2103	08470-2103-01	3-8470210301-7							7.0 lbs
	Size 3	SM 2104	08470-2104-01	3-8470210401-4							8.5 lbs
	Size 4	SM 2105	08470-2105-01	3-8470210501-1							10.10 lbs



Petfine™

Veterinary Products

A line of animal health products, designed to deliver maximum comfort for pets and peace of mind for pet parents.



Product Description	Code	NDC	UPC	Box Qty	Boxes/Case	Item Size	Item Weight	Case Size	Case Weight
Petfine™ Lancets Lancets (28G)	AT 0569VET	08470-0569-01	3-8470056901-3	100	200	2.625" x 3.25" x 1.625"	0.1 lbs	17.25" x 12.75" x 14"	24.3 lbs
Petfine™ Lancing Device Advanced Lancing Device	AT 0275VET	08470-0275-01	3-8470027501-3	1	300	2.125" x 1.375" x 3.625"	0.05 lbs	11.75" x 15.5" x 21"	18.8 lbs
Petfine™ Pentips® 0.25mm (31G) x 8mm Pen Needles with a Built-In Remover	AT 3835VET	08470-3835-01	3-8470383501-6	100	24	4" x 5.25" x 3.625"	0.45 lbs	17.5" x 12" x 12.75"	12.85 lbs
Petfine™ Pentips® 0.33mm (29G) x 12mm Pen Needles with a Built-In Remover	AT 3825VET	08470-3825-01	3-8470382501-7	100	24	4" x 5.25" x 3.625"	0.45 lbs	17.5" x 12" x 12.75"	12.85 lbs
Petfine™ Auto Injector Automatic Syringe Injection Aid	AJ 3310VET	08470-3310-01	3-8470331001-8	1	50	4.875" x 6.75" x 1.5"	0.35 lbs	20.75" x 17" x 10.5"	17.75 lbs

Global Offices



North America

Owen Mumford Inc.
1755 West Oak Commons Court
Marietta, Georgia 30062, USA
T : +1 770 977 2226
F : +1 770 977 2866
E : info@owenmumfordinc.com



Head Office


International Sales & Marketing

Owen Mumford Ltd
Brook Hill, Woodstock
Oxfordshire OX20 1TU, United Kingdom
T : +44 (0)1993 812021
F : +44 (0)1993 813466
E : info@owenmumford.co.uk



UK / Ireland Sales & Marketing

Primsdown Industrial Estate
Worcester Road, Chipping Norton
Oxfordshire OX7 5XP, United Kingdom
T : +44 (0)1608 645555
F : +44 (0)1608 646546
E : info@owenmumford.co.uk



China

Owen Mumford (Shanghai)
Medical Device Company Ltd
25/F Central Plaza, 381 Huai Hai
Middle Road, Huang Pu District
Shanghai, China 200020
T : 021 6032 3528
F : 021 6032 3782
E : info@owenmumford.co.uk



France

Owen Mumford
BP 3333, 27950 Saint Marcel, France
T : +33 (0)2 32 51 88 70
F : +33 (0)2 32 21 96 95
E : contact@owenmumford-fr.com



Germany

Owen Mumford GmbH.
Alte Häge 1
63762 Großostheim-Ringheim, Germany
T : +49 (0)6026/97750
F : +49 (0)6026/977529
E : info@owen-mumford.de



Malaysia

No. 21, Jalan Bioteknologi 3,
Kawasan Perindustrian SILC,
79200 Nusajaya, Johor., Malaysia
T : 60 (0)7 553 9617
E : my.info@owenmumford.com



Mexico

T : 1993 812 021
E : mx.info@owenmumford.com



Making a **World of Difference**
to a World of People

OWEN MUMFORD USA, INC TERMS AND CONDITIONS OF SALE

1. DEFINITIONS

“**Customer**” means the person, firm or company who purchases the Products from OM.

“**OM**” means Owen Mumford USA, Inc.

“**Products**” means any products sold by OM to Customer from time to time.

“**Terms and Conditions**” means the terms and conditions of sale set forth in this document, and no others.

2. APPLICATION OF TERMS

2.1 Subject to Section 2.2, the sale by OM of Products to Customer will be on these Terms and Conditions to the exclusion of all other terms and conditions (including without limitation any terms or conditions which the Customer purports to apply under any purchase order, confirmation of order, specification or other document).

2.2 These Terms and Conditions apply to all OM's sales and supplies of Products unless expressly superseded by a written supply, distribution or other agreement signed by duly authorized representatives of each of OM and Customer, in which case such written agreement shall govern the relationship between OM and Customer as set forth therein. No variation to these Terms and Conditions will be effective unless expressly agreed in writing by a duly authorized representative of OM.

2.3 Each order for Products by the Customer will be deemed to be an offer by the Customer to purchase Products subject to these Terms and Conditions.

2.4 No order placed by the Customer will be deemed to be accepted by OM until OM delivers the Products to the Customer or a written acknowledgment of order is issued by OM. OM reserves the right to charge the Customer for the costs of any materials incurred by OM in relation to a Customer's order following cancellation of that order by the Customer prior to OM's acceptance. Any quotation provided by OM is given on the basis that no contract will come into existence until such acceptance. Any quotation is valid for a period of 30 days only from its date, provided that OM has not previously withdrawn it.

3. DESCRIPTION

3.1 The description of the Products will be as set out in OM's current price list or quotation (as the case may be).

3.2 All drawings, descriptive matter, specifications and advertising issued by or on behalf of OM and any descriptions or illustrations contained in OM's catalogs, brochures or web site are issued or published for the sole purpose of giving an approximate idea of the products described in them. They are not and do not form a part of these Terms and Conditions.

4. DELIVERY

4.1 Unless otherwise agreed in writing by OM, delivery terms of the

Products will be F.O.B. Marietta, GA, to Customers addresses. A ten dollar(\$10.00)fee will be applied to each drop-ship order, in addition to delivery terms of F.O.B. Marietta, GA.

4.2 Any dates specified by OM for delivery of the Products are intended to be an estimate and time for delivery will not be made of the essence by notice. If no dates are so specified, delivery will be within a reasonable time. OM reserves the right to deliver the Products in installments.

4.3 Where the Customer has agreed to purchase a specified quantity of Products, it will give OM reasonable notice of its delivery requirements during that period. Failure by the Customer to take delivery of the total specified quantity of Products in the relevant time period will constitute a material breach hereof.

4.4 If for any reason the Customer will not accept delivery of any of the Products when they are ready for delivery, or OM is unable to deliver the Products on time because the Customer has not provided appropriate instructions, documents, licenses or authorization: (a) all risk in and to the Products, including without limitation all risk of loss, will pass to the Customer; (b) the Products will be deemed to have been delivered to Customer; and (c) OM may store the Products until delivery whereupon the Customer will be liable for all related costs and expenses (including, without limitation, storage, demurrage, insurance and taxes).

4.5 If delivery of the Products is to take place other than at OM's place of business, the Customer will provide at its expense adequate and appropriate equipment and manual labor for unloading the Products.

4.6 If OM delivers to the Customer a quantity of Products more or less than the quantity set out in the relevant order, Customer shall notify OM within two (2) business days of receipt, then: (a) in the case of a shortfall, the Customer shall not be entitled to object to or reject the Products actually delivered by reason of the shortfall; and (b) in the case of a surplus, the Customer shall return such surplus to OM at OM's risk and cost; or (c) issue a purchase order to purchase such surplus Product. No claims for Product shortages will be honored if not communicated to OM within two (2) business days following receipt.

4.7 The Customer is responsible for notifying the courier of any visible damage upon product delivery, noting any discrepancy on the bill of lading. Failure to notify OM and the carrier in a timely manner will result in Products being returned per Section 10.

5. NON-DELIVERY

5.1 The quantity of any Product Shipment as recorded by OM upon dispatch from OM's place of business will be conclusive evidence of the quantity received by the Customer on delivery unless the

Customer can provide conclusive evidence proving the contrary.

5.2 OM will not be liable for any non-delivery of Products (even if caused by OM's negligence) unless written notice is given to OM within 5 business days of the date when the Products would be in the ordinary course of events have been received.

5.3 Any liability of OM for non-delivery of the Products will be limited to, at its option, replacing the Products within a reasonable time or issuing a credit note in respect of such Products.

6. RISK/TITLE

6.1 The Products are at the risk of the Customer from the time of delivery.

6.2 Ownership of the Products will not pass to the Customer until OM has received in full (in cash or cleared funds) all sums due to it in respect of:

(a) the Products; and (b) all other sums which are or which become due to OM from the Customer on any account;

6.3 Until ownership of the Products has passed to the Customer, the Customer will: (a) hold the Products on a fiduciary basis as OM's bailee; (b) not destroy, deface or obscure any identifying mark or packaging on or relating to the Products; (c) maintain the Products in satisfactory condition and keep them insured on OM's behalf for their full price against all risks to the reasonable satisfaction of OM. On request the Customer shall produce the certificate of insurance to OM; and (d) hold the proceeds of the insurance referred to in Section 6.3 (c) on trust for OM and not mix them with any other money, nor pay the proceeds into an overdrawn bank account.

6.4 The Customer may resell the Products before ownership has passed to it solely on the following conditions: (a) any sale shall be effected in the ordinary course of the Customer's business at full market value; and (b) any such sale shall be a sale of OM's property on the Customer's own behalf and the Customer shall deal as principal when making such a sale.

6.5 The Customer's right to possession of the Products shall terminate immediately if: (a) the Customer undergoes any of the events set out in Section 11.2(c); or (b) the Customer encumbers or in any way grants or suffers a lien or security interest on any of the Products.

6.6 OM shall be entitled to recover payment for the Products notwithstanding that ownership of any of the Products has not passed from OM.

7. PRICE

7.1 The price for the Products will be that set out on OM's current price list, acknowledgment of order or quotation (as the case may be). Products are to be purchased in master case quantities, unless otherwise detailed in writing by OM.

OWEN MUMFORD USA, INC TERMS AND CONDITIONS OF SALE

7.2 The price for the Products will be exclusive of any tax.

7.3 OM has a minimum order requirement of one-hundred dollars(\$100.00). Customer orders which fall below the one-hundred dollar(\$100.00) minimum will have a ten dollar(\$10.00) below minimum fee applied.

8. PAYMENT

8.1 OM will invoice the Customer for the Products together with any other amounts pursuant to Section 7.2 upon shipment.

8.2 The Customer will pay the amount invoiced pursuant to Section 8.1 within 30 days of the date of the invoice. Time for payment will be of the essence.

8.3 No deduction whether by way of set-off, counterclaim, discount, abatement or otherwise shall be had by Customer unless otherwise specifically agreed in writing by an authorized representative of OM.

8.4 A finance charge will be leveled at a periodic rate of 1.5% per month (annual percentage 18%) on past due balances over 60 days outstanding.

8.5 Should it be necessary to employ a collection agency or attorney to collect monies due, the Customer will be responsible for any and all reasonable costs of collection.

9. QUALITY

9.1 OM warrants that (subject to the other provisions of these Terms and Conditions) upon delivery, the Products will be suitable for their intended purpose (as set out in OM's catalogs, brochures, or web site (as the case may be)) (the "Warranty").

9.2 OM will not be liable for a breach of the Warranty unless: (a) the Customer gives written notice of the defect to OM, and (if the defect is as a result of visible damage in transit) to the carrier upon delivery; and (b) OM is given a reasonable opportunity after receiving the notice of examining such Products and if requested by OM the Customer returns such Products to OM's place of business at OM's cost for examination.

9.3 OM shall not be liable for a breach of the Warranty if: (a) where the defect is as a result of damage in transit, delivery was at OM's place of business; or (b) the Customer makes any further use of such Products after giving such notice; or (c) the defect arises because the Customer failed to follow OM's oral or written instructions as to the storage, installation, commissioning, use or maintenance of the Products or (if there are none) good trade practice; or (d) the Customer alters, repairs or destroys such Products without the written consent of OM.

9.4 Subject to Sections 9.2 and 9.3, if any of the Products do not conform with the Warranty OM will at its option repair or replace the Products (or the defective part) or refund the price of the Products at the pro rata price therefor provided that, if OM so requests, the

Customer will, at OM's cost, return the Products or the part of the Products which is defective to OM.

9.5 If OM complies with Section 9.4, it will have no further liability for a breach of the Warranty in respect of such Products.

9.6 Any Products replaced will belong to OM and any repaired or replacement Products will be covered by the applicable Warranty.

10. RETURN POLICY

10.1 Prior to returning any Product, Customer must contact OM's customer service for a return authorization, RMA. Any return without an authorization number will be refused. Customer must provide OM with the item number, quantity, and lot numbers of all Product(s) Customer is requesting to return. Return authorization numbers must be indicated on the packing list and written on the outside of all cartons. OM will not accept returns or issue credits for products which are within one (1) year of expiration.

10.2 All Products must be returned by the Customer from which they were originally sold to and be accompanied by their original purchase order number, third party returns will be refused. Upon return to OM, Products must be in good, salable condition, in their original, unopened packaging, with sterility seals fully intact; free from rips, tears, dents, dings, and creases, and showing no signs of being in contact with dust, dirt, debris and/or water.

10.3 Returned goods are subject to OM inspection and determination on condition of the returned Product. If OM determines a Product to not be in good, salable condition, the return may be refused or additional restocking fees may be applied, per the discretion of OM.

10.4 Resalable merchandise must be returned within sixty (60) days of date of invoice. A minimum fifteen percent (15%) restocking fee will apply to all product returns, as determined by OM.

10.5 Return freight charges are the responsibility of returning Customer, unless otherwise detailed in Section 4.7 or 9.4.

11. LIMITATION OF LIABILITY

11.1 Nothing in these Terms and Conditions will operate to exclude or limit the liability of OM for death or personal injury arising out of its negligence or for its fraud in any way whatsoever.

11.2 Subject to Section 10.1, under no circumstances will OM have any liability (whether in contract, tort (including negligence) or otherwise) under or in connection herewith or the sale of Products for: (a) indirect, special or consequential losses of any nature; (b) wasted or lost management time or time of employees, loss of profits, contracts or business, loss of goodwill or loss of anticipated savings; or (c) any increased costs or expenses.

11.3 Subject to Section 10.1, OM's maximum total liability (whether in contract, tort (including negligence) or otherwise) under or in connection with these Terms and Conditions and the sale of the

Products, will not exceed the total price of the Products supplied to Customer.

11.4 The express terms set forth herein are in lieu of all warranties, conditions, terms, representations, undertakings and obligations implied by statute, common law, custom, trade usage, course of dealing or otherwise all of which are excluded to the fullest extent permitted by law.

12. TERMINATION/SUSPENSION

11.1 If Customer fails to make any payment on the applicable due date to OM. OM will be entitled to suspend any further deliveries to Customer without liability to Customer.

12.2 Either party (the "first party") may terminate these Terms and Conditions and related orders with immediate effect by giving written notice to the other (the "second party") without any compensation or damages due to the second party, but without prejudice to any other rights or remedies which the first party may have. If the second party: (a) materially breaches any of these Terms and Conditions which is not capable of remedy; (b) materially breaches any of these Terms and Conditions which is capable of remedy but is not remedied within 30 days of receipt of written notice from the first party specifying the breach and requesting remedy; or (c) has a receiver, administrative receiver, administrator or other similar officer appointed over it or over any part of its undertaking or assets or passes a resolution for winding-up (other than for purpose of a bona fide scheme of solvent amalgamation or reconstruction) or a court of competent jurisdiction makes voluntary arrangement with its creditors or ceases or threatens to cease to carry on business or is unable to pay its debts, or undergoes or is subject to any analogous acts or proceedings under any foreign law, including, but not limited to, bankruptcy proceedings.

13. ASSIGNMENT

Customer may not assign or otherwise transfer any of its rights or obligations hereunder or any part hereof without the prior written consent of OM.

14. DELAY

OM reserves the right to defer the date of delivery or to cancel its agreement with Customer or reduce the volume of the Products ordered by Customer (without liability to Customer) if it is prevented from or delayed in the carrying on of its business due to circumstances beyond the reasonable control of OM.

15. GENERAL

15.1 Customer will comply with all relevant regulatory requirements in connection with the Products including without limitation those relating to licenses, complaints, distribution and storage.

OWEN MUMFORD USA, INC TERMS AND CONDITIONS OF SALE

Further, Customer undertakes to maintain appropriate up-to-date and accurate records to enable the immediate recall of any batches of the Products or any of them from any subsequent purchaser of the Products. These records will include records of deliveries to Customers (including details of batch numbers, delivery date, name and address of customer, and telephone number and fax or telex number if available).

15.2 Each right or remedy of OM hereunder is without prejudice to any other right or remedy of OM whether hereunder or not.

15.3 These Terms and Conditions, together with each relevant offer, quotation, acknowledgment or confirmation constitutes the entire agreement between the parties relating to its subject matter. Each party acknowledges that it has not entered into such agreement on the basis of any warranty, representation, statement, agreement or undertaking except those expressly set out herein.

15.4 If a party fails to enforce or delays in enforcing an obligation of the other party, or fails to exercise or delays in exercising a right hereunder, the failure or delay will not affect their right to enforce that obligation or constitute a waiver of that provision on a future occasion.

15.5 No person who is not a party hereto has any right to prevent the variation or cancellation of any provision hereof or its termination, and no person who is not a party hereto may enforce any benefit conferred upon them hereby, unless expressly provided otherwise.

15.6 These Terms and Conditions will be governed by and construed in accordance with the laws of the State of Georgia, without application of its conflicts of laws rules.

16. NOTICES

16.1 Any notice between the parties hereunder must be in writing and delivered by hand or sent by pre-paid first class post, or courier facsimile transmission; or (a) (in case of communications to Customer) to the registered office of the address as shall be notified to Customer by OM; or

(b) (in case of communications to Customer) to the registered office of the addressee (if it is a company) or (in any other case) to any address of Customer set out in any communication between Customer and OM or such other address as shall be notified to OM by Customer.

16.2 Communications shall be deemed to have been received: (a) if sent by pre-paid first class post, 5 business days; (b) if delivered by hand or courier on the day of delivery; (c) if sent by facsimile transmission on a working day prior to 4:00pm, at the time of transmission and otherwise on the next working day.